



UTILITY PATENT APPLICATION TRANSMITTAL Submit an original and a duplicate for fee processing (Only for new popprovisional applications under 37 CFR 1.53(b))

	(Only for new horiprovisional applications under 37 CFR 1.53(b))					
ADDRESS TO: Assistant Commissioner for Patents Box Patent Application Washington, D.C. 20231		Attorney Docket No. MBHB00-203 First Named Inventor Ruderman, Frank R. Express Mail No. EL442910428US Total Pages				
	APPLICATION ELEMENTS	ACCOMPANYING APPLICATION PARTS				
1. 2. 3. 4. 5.	Transmittal Form with Fee Specification (including claims and abstract) [Total Pages 25] Drawings [Total Sheets 25] Oath or Declaration[Total Pages 2] a. unexecuted b. Copy from prior application [Note Boxes 5 and 17 below] i. Deletion of Inventor(s) Signed statement attached deleting inventor(s) named in the prior application Incorporation by Reference: The entire disclosure of the prior application, from which a copy of the oath or declaration is supplied under Box 4b, is considered as being part of the disclosure of the accompanying application and is hereby incorporated by reference therein. Microfiche Computer Program Nucleotide and/or Amino Acid Sequence Submission a. Computer Readable Copy b. Paper Copy c. Statement verifying above copies	8. Assignment Papers 9. Power of Attorney 10. English Translation Document (if applicable) 11. Information Disclosure Statement (IDS)				
17.	If a CONTINUING APPLICATION, check app	propriate box and supply the requisite information:				
☐ Continuation ☐ Divisional ☐ Continuation-in-part of prior application Serial No.						

	5 (2m);	AP	PLICATION FEES		
BASIC FEE	ي ميان در	· 12		30	\$690.00
CLAIMS	NUMBER F	ILED	NUMBER EXTRA	RATE	
Total Claims	20	-20=		x \$18.00	\$
Independent Claims	2	- 3=		x \$78.00	\$
☐ Multiple Dependent Claims(s) if applicable +\$270.00					\$
Total of above calculations =					\$
Reduction by 50% for filing by small entity =					
Assignment fee if applicable + \$40.00					\$
# 55 C	2,0302 2,59°			TOTAL =	\$ 690.00



	TAPPLICATION TRANSMITTAL	Attorney Docket No. MBHB00-203			
18. Please	18. Please charge my Deposit Account No. 13-2490 in the amount of \$.				
19.	19. A check in the amount of \$ is enclosed.				
 20. The Commissioner is hereby authorized to credit overpayments or charge any additional fees of the following types to Deposit Account No. 13-2490: a. ☐ Fees required under 37 CFR 1.16. b. ☐ Fees required under 37 CFR 1.17. c. ☐ Fees required under 37 CFR 1.18. 					
21. The Commissioner is hereby generally authorized under 37 CFR 1.136(a)(3) to treat any future reply in this or any related application filed pursuant to 37 CFR 1.53 requiring an extension of time as incorporating a request therefor, and the Commissioner is hereby specifically authorized to charge Deposit Account No. 13-2490 for any fee that may be due in connection with such a request for an extension of time.					
	22. CORRESPONDENCE ADD	RESS			
Name McDonnell Boehnen Hulbert & Berghoff					
Address	32 nd Floor, 300 South Wacker Drive				
City, State, Zip	Chicago, Illinois 60606				
23. SIGNATURE OF APPLICANT, ATTORNEY, OR AGENT REQUIRED					
Name	Robert J. Irvine III				
Signature	Polut Di I				
Date	March 24, 2000				

UTILITY (Rev. 11/18/97)

CERTIFICATE OF MAILING BY "EXPRESS MAIL" (NEW PATENT APPLICATION)

Express Mail No. EL442910428

Deposited March 24, 2000

I hereby certify that the attached correspondence, identified below, is being deposited with the United States Postal Service as "Express Mail Post Office to Addressee" under 37 CFR § 1.10 on the date indicated above and is addressed to the Asst. Commissioner for Patents, Box New Patent Application, Washington, DC 20231.

By: (person actually depositing)

Patent Application of: Ruderman, Frank R. et al.				
Title: Cardiovascular Healthcare Management System and Method				
X Patent Application (25 pages, including claims)				
X Drawings (25 sheets)				
X Postcard				
X Transmittal Letter				
Check				
Information Disclosure Statement				
Small Entity Form				
Preliminary Amendment				
X Patent Express Mail Certificate				

Attorney Docket No.: MBHB00-203

CARDIOVASCULAR HEALTHCARE MANAGEMENT SYSTEM AND METHOD CASE No. MBHB00-203

Inventor:

Frank R. Ruderman

248 Bayview Drive

San Carlos, California 94070

David T. Shewmake, Ph.D.

419 Baden Street

San Francisco, California 94131-2831

Assignee:

Berkeley HeartLab, Inc.

1875 South Grant St., Suite #700

San Mateo, CA 94402

5

BACKGROUND OF THE INVENTION

This application claims the benefit under 35 U.S.C. § 119(e) of U.S. Provisional Patent Application Serial Number 60/126,003, filed March 24, 1999, entitled "Method for Individualized Patient Treatment" for all common subject matter disclosed therein, and the contents of said application are hereby incorporated by reference. This application also claims the benefit under 35 U.S.C. § 119(e) of U.S. Provisional Patent Application Serial Number 60/168,354, filed December 1, 1999, entitled "Method for Individualized Patient Treatment" for all common subject matter disclosed therein, and the contents of said application are hereby incorporated by reference.

10 1. <u>FIELD OF THE INVENTION</u>

This invention is in the field of managing patient health through an automated online database. More particularly, the system provides for diagnosis and management of cardiovascular disease and includes assessment of cardiovascular risk factors and providing patient treatment plans.

15 2. <u>BACKGROUND OF THE ART</u>

The art describes cardiovascular risk factors such as age, smoking, weight, family history, blood pressure, lipid profiles including low density lipoprotein (LDL) and high density lipoprotein (HDL) and subclasses (fractions) of LDL and HDL. Methods for measuring these factors and relating them to patient treatment are also known. Generally, physicians assess a patient's risk factors, make a diagnosis based on test results and symptoms and manage patient treatment through drugs, exercise, diet and a variety of surgical techniques. Instructions are generally given directly to the patient by the physician. Patient compliance generally involves interview in follow-up office visits.

SUMMARY OF THE INVENTION

The present invention provides a cardiovascular healthcare management system that has an infomediary site, which includes databases with cardiovascular disease risk assessment and cardiovascular disease management information. Patient data is stored in an information database, and patient test results are stored in a clinical database, both of which are selectively accessible to both physicians and patients. Patients' test results and

personal information are thus added to the databases and may be viewed by the physician to assist in diagnosis.

A knowledge base may also be included to provide further diagnostic data to a physician. The knowledge base may be used as part of a diagnostic engine that analyzes test results to determine whether they meet certain criteria indicative of a patient's health condition, or might simply contain descriptions of risk factors or combinations of risk factors that are believed to be indicative of specific health conditions. The knowledge base of the preferred embodiment includes software referred to herein as a diagnostic engine that performs risk assessments on specific test results and stores risk assessment information in the database for later viewing by the physician to further assist in diagnosis.

The database also has information related to managing the patient's cardiovascular health. The cardiovascular disease (CVD) management system provides for the physician having electronic access to the infomediary site for receiving patient test results and entering diagnosis and treatment information to facilitate the building of a treatment plan. The treatment plan is created after viewing test results and relevant patient data and entering any supplemental information including a diagnosis. The treatment plan may include a recommended diet, prescription (and nonprescription) drugs, an exercise regimen, and alternate cardiovascular products that may be available for purchase through the infomediary site (e.g., blood pressure cuffs to monitor blood pressure, anti-embolism support hosiery, dietary products, educational materials, etc.). Portions of the treatment plan are preferably derived from templates that are provided by the infomediary site 100. In addition, the infomediary site may present templates having the most relevant treatment plan components based in part upon the physician's previously entered diagnosis information.

The system also provides for electronic communication between the patient and physician by way of access to one or more records in the infomediary site. In this way, a physician can electronically provide cardiovascular disease management instruction such as diet, exercise and medication to the patient and the patient can electronically provide compliance information to the physician directly or indirectly via an alternate input device.

20

5

In another aspect of the infomediary site, a case manager may be provided to collect and input initial patient data and to perform follow-up services like collecting and inputting compliance data to the infomediary site. The case manager may also review initial patient history, test results, and then create a suggested treatment plan for the cardiologist to review. The cardiologist may accept the treatment plan unchanged or modify it, and then electronically release the treatment plan to the patient's folder, which stores personalized information within the CVD management system site. The case manager also proactively contacts patients to advise and assist in the implementation of their treatment plan.

Alternatively, tele-medicine devices may be used to gather compliance information for submission to the infomediary site.

BRIEF DESCRIPTION OF THE FIGURES

Figure 1 is a block diagram of the infomediary site;

Figure 2 is a login screen;

Figure 3 is the physician's welcome screen;

Figures 4-12 show additional screens available to the physicians;

Figure 13 is the patient's welcome screen;

Figures 14-17 show additional screens available to the patient;

Figures 18-22 show the database tables of the preferred embodiment;

Figure 23 is a schematic of the cardiovascular informatics;

Figure 24 is a more detailed block diagram of the infomediary site; and

Figures 25 and 26 are flow diagrams for methods associated with the infomediary site.

DETAILED DESCRIPTION OF THE INVENTION

Figure 1 illustrates a schematic of the infomediary site <u>100</u>. It shows physician <u>102</u> in communication with the site <u>100</u> to retrieve test results from a database <u>104</u>, communication of diagnosis information and a treatment plan or other data to the database <u>104</u> by way of server <u>112</u>. The server scripts <u>112</u> generally provide a "front end" to access the records within database <u>104</u>. Specified records of database <u>104</u> are also made accessible to the respective patients <u>106</u>. The patient <u>106</u> can view test results,

15

20

25

30

treatment plans, and communications from the physician <u>102</u>, and communicate compliance information to the physician <u>102</u>. The lab <u>108</u> provides for inputting test results into the site. There may also be an optional guest <u>110</u> access for reviewing limited data in the system.

5 PHYSICIAN ACCESS

The physicians 102 can log-on to the infomediary site 100 and communicate with patients 106 and the infomediary site 100 to manage the patient's cardiovascular healthcare. The physician 102 first logs in to the site 100 at the login screen 150 by entering a username in field 152, a password 154, and submitting the data by clicking on submit button 156. The physician is then presented with the welcome screen 160 as shown in Figure 3. The interface made up of the web pages that provide physician access to the site is referred to herein as the physician data access interface. The physician welcome screen 160 includes a navigation bar having a plurality of links represented by buttons 161-169. The physician welcome screen also includes summary views of certain data such as the physician's statistics at tab 174, patient activity at tab 176 and professional news at tab 178. Tabs 174, 176, and 178 are also links to more detailed screens, which can also be accessed by buttons 164, 161, and 168, respectively. Preferably, the physician is able to customize his or her own welcome screen to provide customized summary information.

The link underlying patient activity button 161 brings the physician to the patient activity screen 200 as shown in Figure 4. The patient activity screen 200 is generated by ASP script physician_patient_activity.asp. The patient activity screen 200 provides a list of patients who may need attention. The physician's attention may be needed due to the receipt of new test results (input by the Lab 108), the receipt of patient correspondence or compliance data, or the like. The link underlying patient list button 162 brings the physician to the patient list screen 210 as shown in Figure 5, which is preferably generated by the ASP script physician_patient_list.asp. The screen 210 contains a list of all current (and/or past) patients for that particular physician. The link underlying e-mail button 163 brings the physician to the email screen 220 as shown in Figure 6, which is preferably generated by ASP script physician_email_detail.asp.

The names of one or more of the patients of the particular physician are displayed on the screens in Figures 4-6. Preferably, the patient name (e.g., 201, 211, and 221) is also a hypertext link to the patient's data file, thereby providing the physician with convenient access to the detailed patient record displayed on a patient folder screen, one of which is shown at 230 in Figure 7. The patient folder 230 is preferably generated by ASP script physician patient activity details.asp.

The patient folder <u>230</u> has a general patient information display area <u>232</u>, patient correspondence display area <u>234</u>, test results display area <u>236</u>, FAQ audit trail display area <u>238</u>, and compliance data display area <u>240</u>. The patient folder may also contain patient billing and insurance coverage and/or insurance carrier information (not shown). The physician can create initial patient histories by utilizing and editing standard forms. Alternatively, a patient may provide medical background information to the infomediary site <u>100</u>, either by way of the network through an initial registration procedure, or by calling a data entry specialist provided by the infomediary site <u>100</u>. Preferably, the infomediary site provides the physician with the services of the case manager, who is typically a lipid nurse, nurse, physician assistant, or other qualified person. The case manager may also act as the data entry specialist.

The general patient information display area <u>232</u> contains general patient information including name, birth date, medications, insurance information, and salient medical history information, etc.

The patient correspondence display area <u>234</u> is referred to herein as a communication interface. It is presented as an email inbox, but preferably does not operate as standard internet email. That is, physicians preferably communicate with their patients by way of the database <u>104</u> located at the infomediary site <u>100</u>. Referring back to Figure 1, the physicians preferably enter data to be communicated into one or more records of database <u>104</u> by way of the infomediary server scripts <u>112</u>. The server scripts <u>112</u> process physician and patient "email" entries to the database <u>104</u> and provide an email-type interface. Preferably, actual emails are not sent by the physician, patient or infomediary site <u>100</u>. Standard email may be utilized in an alternative embodiment; however, less monitoring and control (e.g., delivery verification) of the communication is available when standard email is used. To generate an internet email from an ASP script,

10

15

20

25

30

for example, an ASP server running with the NT 4.0 Option Pack installed and Microsoft SMTP Service (from the Option Pack) running may use the CDONTS component. The correspondence display area 234 also shows the status of the email with icons 235, as well as an indication of whether the physician has replied to the patient's email. Of course, the infomediary site may utilize standard internet email to notify the patient of upcoming appointments for lab test work, doctor office visits, compliance conferences with case managers, or to request that the patient log in to the infomediary site to view newly posted treatment plans, compliance requests, or for other similar notifications.

The FAQ audit trail display area 238 allows the physician to view which general information the particular patient has accessed over the network. This may reveal additional information to the physician about the patient's concerns that the patient would not otherwise convey to his physician. The compliance data display area 240 provides the physician with information provided by the patient as to whether the patient has been complying with his or her treatment plan. Preferably, the patient provides information to the case manager who, in accordance with the treatment plan, initiates follow-up contact with the patient. The case manager solicits information regarding the patient's compliance with the treatment plan such as diet, exercise, medication, and weight and blood pressure, if available. The case manager may input the compliance data into the information database 104. Alternatively, the patient may input the data directly via the infomediary site <u>100</u>. In addition, with reference to Figure 24, the patient may utilize tele-medicine devices to provide compliance-type data. For example, health related measurement devices such as a weight scale 115, blood pressure cuff 116, electrocardiogram devices 117, and the like, may be connected to the network 101 for transmitting data to the infomediary site 100. The device may be connected directly to the network <u>101</u>, or via the patient's home computer, or through a networking device in the patient's home. The tele-medicine devices may also be made available at another location such as a drug store, etc. The device preferably communicates the measured data and automatically submits the data or compliance information to the infomediary site. The tele-medicine device may include software (i.e., a resident device driver, or java-type plug in module) to format appropriately the identification and authentication information along with the data for transmission to the site 100. The physician may be

10

15

20

25

30

notified of the presence of compliance data either via the email inbox $\underline{234}$ or by listing that patient in the patient activity display screen $\underline{200}$.

The test results display area 236 shows recent test work, and a test results link 237 to allow the physician to view the test results. A typical test results screen 260 is shown in Figure 8. In particular, the test results of Figure 8 show the patient's HDL Gradient Gel Electrophoresis in tabular form by HDL subclass, and in graphical form. As seen in the table, the test results show the normalized percentages of HDL subfractions, as well some points (10%, 50% and 90%) on the distribution of normalized percentages. Specifically, the table shows that 10% of people in the sample have a normalized percentage of HDL Subclass 2b of 4% or more, while 90% have a normalized percentage of 18% or less, with the 50% point being 8%. Thus, for example, the physician can see that the particular 2b Subclass test result of 15% is generally within the higher ranges. It should be noted that more complete tables may be provided to show finer resolution of percentages. It should also be noted that the people in the sample to which the test results are compared may vary. In particular, it is often helpful for a physician to be able to compare test results to a population having a similar diagnosis. To illustrate this, in the example shown in Figure 8, the percentages listed are for those people diagnosed with coronary artery disease (CAD). Thus, the physician is able to determine that the particular results show levels of HDL Subclass 2a and 2b that are higher than average compared to people with CAD, and levels of HDL Subclass 3a, 3b and 3c that are lower than average compared to that same population, and interpret the results accordingly. In a preferred embodiment the physician may dynamically select the reference population to further assist in making a diagnosis. The doctor preferably has the ability to adjust a number of analysis parameters to customize the test result analysis or comparison. The parameters may specify aspects of the general population used as a reference and may include age range, diagnoses, symptoms, LDL and HDL subfraction ranges, etc. The site preferably provides default ranges of many analysis parameters to yield comparisons that have been determined (via research, data mining, etc.) to be useful in the diagnosis or treatment of CVD.

The physician can also view and print current and prior test results including patient personal information - name, address and telephone number. With the

10

15

20

25

30

infomediary site <u>100</u>, the physician can view cardiovascular risks from the knowledge base, including graphic representation of the patient's condition and test results. Charts or histograms comparing current and prior test results are also available from the infomediary site database.

The patient folder <u>230</u> also includes a diagnosis link <u>250</u>, treatment plan link <u>252</u>, and letter of consultation link <u>254</u>. The diagnosis link <u>250</u> retrieves the diagnosis screen <u>290</u> as shown in Figure 9. The various diagnoses are shown, with the dates and summary information. For example, Diagnosis 1, made in 8/99, would be displayed in field <u>291</u>. The physician may also create a new diagnosis using link <u>292</u>.

The treatment plan link 252 retrieves treatment plan screen 300 as shown in Figure 10. The physician can create treatment plans for their patients from information in databases of the infomediary site <u>100</u>. Treatment plan screen <u>320</u> includes overview display area 302, exercise display area 304, prescription medication display area 306 (generic and brand names, as well as the ability to print the prescription, or to specify a pharmacy to where the prescription may be electronically forwarded), diet display area 308, cardiovascular (CV) products display area 310, suggested readings display area 312, adjunct medications display area <u>314</u>, and follow up display area <u>316</u>. Each of the display areas 302 – 316 include links 320 to pages having further details relating to the respective sections. Preferably, the infomediary site provides the physician with a treatment plan template having the above-referenced components. The physician may select specific entries by a number of methods, including drop-down boxes, text boxes, bullet lists, etc. In addition, the template may initially provide selections that are consistent with the patient's data, including one or more of the following: test results, most recent test results, salient patient history, physician diagnosis data, etc. The finalized treatment plan may also be reviewed by the infomediary site to check for consistency of the plan with the same factors, and the physician's attention may be directed to aspects of the treatment plan that appear to conflict with the particular patient's data.

The link <u>254</u> to the letter of consultation allows the physician to construct letters of consultation from standard forms, patient data, standard paragraphs and phrases in the

10

15

20

25

30

infomediary site database <u>104</u>. These letters may be printed, stored and sent electronically or through regular mail.

The physician can also, through the infomediary site, visit e-commerce vendor's websites to purchase and recommend purchasing to their patients of prescription medication, adjunct medication, exercise equipment, dietary products, cardiovascular products and educational materials, such as reading materials, software, video recordings, etc.

The infomediary site 100 also provides physicians with educational materials such as recent news and suggested readings. The infomediary site 100 also provides physicians with administrative function capability. The physician can perform administrative functions, including creating or editing a physician profile, managing billing information and logs, and viewing daily activity summaries.

The infomediary site 100 may also perform insurance pre-certification to ensure that requested tests are covered by the patient's insurance. If the test is not covered, the physician and/or patient is notified, and further instructions are requested as to whether the physician and/or patient nevertheless wishes to have the test performed. The communication can be via telephone from the case manager, or via the infomediary site through the physician's "email" inbox. Typically, the infomediary will telephone the patient to inform them that the test is not covered or partially not covered. The nature of the test and the reason for the test is explained to the patient, and the patient may then elect to pay for the scheduled test or to decline the test. In the event the test is declined, the requesting physician is notified.

PATIENT COMMUNICATION

Patients can log-on to the infomediary site 100 and view screens that provide links to other pages or news. The patient 106 first logs in to the site 100 at the login screen 150 (Figure 2) by entering a username in field 152, a password 154, and submitting the data by clicking on submit button 156. The patient is then presented with the patient welcome screen 340 as shown in Figure 13. The patient welcome screen 340 includes a navigation bar having a plurality of links represented by buttons 351-359. The patient welcome screen also includes summary views of certain data such as test results 341 (as shown further in Figures 14 and 15), diagnosis information 342, treatment

10

15

20

25

30

plans <u>343</u> (as shown further in Figure 16), and email inbox <u>344</u>. Figure 17 depicts a compliance input screen that may be accessed directly by the patient, or by the case manager.

The patient can communicate with their physician through a communication interface. Preferably, this takes the form of mutual access to records in the infomediary site to obtain doctor treatment plans and report their compliance to the doctor as described above with respect to physician email. The patient can purchase products, including medications, adjunct medications, exercise equipment, dietary products, cardiovascular products and educational material directly from the infomediary site <u>100</u>. Additionally, or alternatively, the infomediary site may contain links to vendor websites for purchasing the products from a third party.

NON-DOCTOR/NON-PATIENT ACCESS TO THE INFOMEDIARY SYSTEM

The infomediary system is adapted so that patient results from laboratory tests and related information can be inputted into the system. In particular, the site 100 may be accessed by the lab 108 so that test result data may be input directly through the web site server scripts. Alternatively, test result data may be input to the system by a data entry specialist or by a case manager.

The infomediary site may also be accessed by guests to the system. Guests may browse information, FAQ lists, physician profiles, physician referral services, or request further information or follow-up contact.

CARDIOVASCULAR INFORMATICS

The infomediary site has advanced cardiovascular informatics related to factors such as genetics, diet, exercising and medication. The cardiovascular informatics is illustrated in Figures 23 and 24. Table 902 illustrates lipid-related data in the infomediary site 100. The relationship of lipids 902 to genetic markers 900 and other coronary risk and disease factors such as those listed in 904 are analyzed by and stored in the CV Therapeutic database 906. Of course, as an initial matter, an individual patient might not have every test listed in 900 and 902; however, the infomediary site includes a diagnostic engine 113 to identify whether any additional tests should be recommended to a physician. Note that tables 900, 902, 904 and 908 are shown in summary form, but

10

15

20

25

30

generally take the form of relational database tables such as those shown in Figures 18-22, having ID fields, text fields, patient and test ID fields, etc.

Thus, a cardiovascular risk database is formed from the information in tables 900, 902, and 904 within clinical database 105. Each entry in tables 900 and 902 also includes fields to identify the patient ID, physician/client ID as necessary, and may optionally include one or more suggested diagnosis fields (not shown). The diagnosis fields may also be included in another table associated with the patient within clinical database 105 or information database 104. The diagnosis fields preferably contain diagnosis ID's corresponding to the diagnosis text shown in table 904. The diagnosis fields are populated by the diagnostic engine 113, and are used to convey the analytical results of the diagnostic engine to the physician through the physician access interface. The CV therapeutic database 906 includes algorithms to identify patterns of risk factors, and may compare inputted patient results for lipid subfractions to the normal values. The algorithms are shown in Figure 24 as a diagnostic engine 113. These algorithms may take the form of software programs or scripts such as SQL scripts or the like. The diagnosis engine 113 preferably operates on new test results to identify those results (in combination with other risk factors) that fall within the parameters specified by the algorithm. Additionally, the CV therapeutic database 906 may analyze the lipid subfractions in combination with other patient data such as that listed in table 900 and 904. The results, including possible treatments 908, can then be related to the physician.

As described above, the measurement of the patient's lipids, genetic markers and disease and disorder risk factors and comparison to the database provides a patient's specific risk factors which can be related to possible treatments. Furthermore, the book "The Heart Disease Breakthrough" by Thomas Yannios, M.D. John Wiley & Sons in 1999 illustrates additional important principals in cardiovascular healthcare that are incorporated into the CV therapeutic database 906 and diagnostic engine 113. This book is incorporated herein by reference. The book points out the important role of LDL and HDL subfractions in cardiovascular disease. This book further discusses the role of genetics, diet, exercising and medication in cardiovascular disease management.

A doctor reviewing the patient's risk factors can formulate an individual course of treatment, and the infomediary site 100 provides a list of possible treatments 908,

10

including exercise, diet, pharmaceuticals, which are preferably customized by the CV Therapeutic database 906 in response to the risk factors identified by the test results of the patient. The physician sends an individualized treatment plan to a patient record in the infomediary site and the patient sends compliance information to the physician to the patient file. Either may access the relevant treatment plan records, but preferably only the physician may modify the recommended plan, and the patient may add, update or modify the compliance data.

The quantitative measurement of LDL and HDL lipid fractions and other lipid measurements are important patient data that is evaluated in the infomediary site. U.S. Patent 5,925,229 describes methods for quantitating LDL fractions for evaluating cardiac disease risk and that patent is incorporated herein by reference. Typical lipids and their method for measurement are listed in Table I.

TABLE I

	<u>Analysis</u>	Methodology	Summary Description of Assay		
5	Apoprotein A1	Immunoturbidimetric Analysis	An insoluble turbid immunoprecipitate is formed by the reaction between Apo A-1 antigen in human plasma and specific antibodies. Resulting turbidity is measured spectrophorometrically. Concentration of Apo A-1 is determined from a five-point calibration curve. Patients run in duplicate.		
10	Apoprotein B	Immunoturbidimetric Analysis	An insoluble turbid immunoprecipitate is formed by the reaction between Apo B antigen in human plasma and specific antibodies. Resulting turbidity is measured spectrophotometrically. Concentration of Apo B is determined from a five point calibration		
15			curve. Patients run in duplicate.		
	Apoprotein-E Isoforms	Isoelectric Focusing	Isoelectric focusing of plasma sample is followed histochemical visualization of the Apo-E bands.		
20	Lipoprotein (a)	ELISA= Enzyme Linked Immunoobsorbant Assay	Monoclonal/Polyclonal sandwich assay. Monoclonal antibody coated wells are used to capture Lp(a) from the sample. Polyclonal anti-Lp(a) horseradish peroxidase conjugate reacts with substrate and a		
25			chromogen producing colored solution. Concentration of Lp(a) mass (mg/dl) quantitatively determined by standard curve		
30	Lipid Profile Cholesterol Triglyceride HDL Cholesterol	Enzymatic Enzymatic Dextran & Magnesium sulfate	Patients run in duplicate: Photometric quantification Photometric quantification Precipitation of LDL and VLDL from plasma followed by assay of		
35	LDL Cholesterol	Calculated - Friedewald Formula	supernatant for cholesterol. LDL=T Chol – HDL Chol – VLDL (Trig/5)		
	Homocysteine	FPIA = Fluorescence Polarization Immunoassay	Specific monoclonal antibody detection with fluorescently labeled analog tracer detection.		
40	LDL GGE Gradient Gel Electrophor		Plasma low density lipoproteins (LDL) are characterized for particle size and distribution of lipid stained mass using polyacrylamide gradient gels (2-14%). LDL particle size analysis is based on a four-point curve using lipoprotein standards correlated by		
45			analytical ultracentrifugation (ANUC). Results of phyenotype analysis are validated by comparison of GGE analysis with ANUC.		
50	HDL GGE	Gradient Gel Electrophoresis	Plasma high density lipoproteins (HDL) are characterized for particle size and distribution of protein stained mass using polyacrylamide gradient gels (3-31%). HDL particle size analysis is based on calibration with globular protein standards. Results validated by comparison of GGE analysis with		
55			analytical ultracentrifugation (ANUC).		

10

15

20

25

30

Patient data from lipid analysis and other data have been utilized to develop a knowledge database within the infomediary site. This knowledge base is useful in analyzing an individual patient's data and also identifying new and previously unknown relationships between test results that aid in the diagnosis and treatment of patients.

Additional relationships are identified from the tables <u>900</u>, <u>902</u>, and <u>904</u> using standard data mining techniques, including genetic algorithms, decision tree induction, association discovery, fuzzy logic, etc. The relationships are then incorporated into the CV therapeutic database <u>906</u>. More particularly, the relationships may be incorporated into programs or scripts (e.g., SQL scripts) within the diagnostic engine <u>113</u>.

Considering 954 patient samples (458 cases and 496 controls), age was a very significant predictor of CVD. Cases are significantly older than controls (60 vs. 52 years of age). After adjusting for the age difference, none of the risk factors are significantly different between the cases and controls. Thus, using the 954 patients, all of the differences in risk factors that exist between cases and controls are due to age, not disease status. All of these patients are high risk and the younger patients have not yet shown clinical manifestation of cardiovascular disease.

Using a subset of age matched cases (N=173, means 60 yr.) and controls (N=173, means 59 yr.) between 54 and 66 years of age, the <u>cases</u> had significantly:

Higher homocysteine (9.7 vs. 8.7, P<0.01), and

Lower TC (179 vs. 201, p<0.0001), LDLC (107 vs. 121, p<0.001), triglyceride (140 vs. 163, p<0.05), apoA1 (112 vs. 123, p<0.01) apoB (85 vs. 96, p<0.001), and TC/HDL2b (14.8 vs. 20.2, p<0.05). These data indicate that the cases are more aggressively treated with medications than the controls.

Using a <u>subset of age-matched</u> cases (N=146, mean 55 yr.) and controls (N=93, mean 55 yr.) between 44 and 66 years of age <u>without hyperlipidemia</u>, the <u>cases</u> had:

Higher HDL3b (19.9 vs. 17.9, p<0.05), HDL3 (58.8 vs. 55.7, p=0.08) and LDLII+IV/HDL2+3 (0.40 vs. 0.38, p=0.11), and

Lower TC (182 vs. 205, p<0.001), LDLC (109 vs. 124, p<0.01), HDLC (44 vs. LDL11A (16.8 vs. 18.2, p=0.09), HDL2b (15.5 vs. 18.6, p<0.05), and HDL2 (41.3 vs. 44.5, p=0.06). These data again indicate that cases may be more aggressively treated with medications than the controls, even though they do not have hyperlipidemia. These

10

15

20

25

30

data also indicate some important risk factors in the cases: a higher ratio of small LDL to HDL, small LDL size and lower HDL2b.

These data illustrate the value of the cardiovascular informatic knowledge base in deriving heretofore unrecognized relationships between data, especially highly discriminating lipoprotein subfractions, in diagnosing risk factors which may govern the treatment of patients.

DETAILED DESCRIPTION OF THE INFOMEDIARY SITE

In a preferred embodiment, the infomediary (information intermediary) site 100 is accessed over a network 101, either a public network such as the internet, or a private network such as a private LAN/WAN, either of which may be accessed over high speed digital subscriber lines or optional dial up (circuit switched) connections provided by a remote authentication server or the like. The network 101, as shown in Figure 24, may include landline and/or wireless access (typically CDMA), as is known in the telecommunications art.

The healthcare information is typically displayed on a web browser running on a personal computer or other web-content viewing device. In the case of patient access, this may include a personal computer 107 with internet access, a television set with a settop decoder (such as that provided by WorldGate Communications, Microsoft's WebTV, OpenTV interactive television, or the like), and may also include a cable modem or wireless modem, or may be a smaller internet appliance device. The patient may also utilize tele-medicine devices such as a weight scale 115, blood pressure cuff 116, electrocardiogram device 117 or other similar monitoring devices.

In the case of a physician, any of the above access devices may be used. Preferably, physicians utilize a portable device such as a clipboard or tablet computer 107 such as those available from Fujitsu, a Clio device available from Vadem, or a smaller device such as a Palm Pilot. The device may have wireless access 111 to the network for data transfers, or may store information locally on a computer 103 and periodically synchronize with the infomediary site via the network.

The information pages are conveyed to the browsing device over the network in the form of a device independent language such as hypertext markup language (html) using hypertext transfer protocol (http). These protocols are typically carried by a lower

level transfer protocol such as those associated with the TCP/IP protocol suite. The device independent language is preferably generated dynamically by the infomediary site, using Common Gateway Interface (CGI) scripts, Active Server Pages (ASP), or other content generation languages and protocols. CGI and ASP are language-independent frameworks for coding of server-side scripts that are executed by Web server 112 in response to a user's request for a universal resource locator (URL). The infomediary site may also be accessible from an interactive television channel having regular multimedia displays of health-related educational information, documentaries, and the like, with the healthcare management features described herein being accessible to patients using software or scripts written in accordance with known interactive television transmission platforms and protocols, such as that provided by OpenTV's suite of EN2 set-top box software, software development tools, studio authoring tools, and Openstreamer data streaming software.

In response to user requests, the infomediary server scripts <u>112</u> access the database <u>104</u>, which contains numerous data tables as shown in Figures 18-22. For example, ASP scripts provide database access via Microsoft's ActiveX Data Objects (ADO), which allows access to ODBC or OLEDB compliant data source including Microsoft Access (Jet), Microsoft SQL Server, and Oracle databases.

The data tables will now be described with reference to Figures 18-22. The various tables in Figures 18-22 make up a relational database, where the tables are related as indicated by the interconnecting lines. As is known to those of skill in the art, tables within relational databases may exhibit a one-to-many relationship. An exemplary description of such a relationship will be given with reference to the UserID of the TWebUsers table 800. The remaining tables are related in a similar manner as is apparent from the drawings.

Table 800 contains a record for each of the users of the infomediary site 100. Each UserID is a unique number defining a particular user's characteristics, including their login name, password, usertype, and a timestamp indicating when the record was created. The tables 801, 803, and 805, all include fields that are linked to the UserID within table 800. The key symbol at the end of the lines connected to table 800 indicate that the value within the UserID field specifies a unique record within the tWebUsers

table <u>800</u>, while the infinity symbol on the other end of the lines connected to tables <u>801</u>, <u>803</u>, and <u>805</u> indicate that more than one record in those tables can be linked to the same UserID in table <u>800</u>. As is know to those of skill in the art, relating the tables in this manner results in a more efficient and powerful database structure than a single large data table.

The table contents will now be generally described, in view of the above-defined relationships. When a patient initially registers with the infomediary site 100, he or she is associated with a physician. The registration is typically performed by the physician's office assistants (independently or with the aid of infomediary site personnel) in advance of the patient's first visit to the infomediary site. With reference to Figure 20, client patient table 840 links a patient's webuser ID to the relevant physician webuser IDs.

Table 841 links physicians (by the clientID) to a physician profile including an affiliated hospital, degree, and region (further described in client region table 844). As shown in client assistant map table 842, each physician assistant has an assistant ID, a clientID (i.e., physician), and a permision level (as described in permission level table 843).

Physicians may be categorized according to client category map table 846 (as described by client news category table 847, and news stories may be categorized as shown in table 845.

With reference to Figure 18, tables 801, 802, 803, 804, 807 and 808 are used to gather the patient history data. Either the patient or data entry personnel (such as a physician's assistant or a case manager or technician provided by the infomediary site) using an appropriate interface (preferably a GUI, such as a web browser, or another database entry interface) collects data that is inserted into the tables. The tables are preferably accessed using a data entry interface such as the web pages provided by the web server 112. When the patient or data entry personnel navigate to the relevant input screens, the database management system 115 (or server scripts 112, such as ASP server scripts) queries the question map table 801 to determine which questions the physician wishes to ask his patients. Each physician (identified by the ClientID field in table 801, which is in turn linked to a UserID in table 800) may create a customized list of questions from a question pool. The list of questions is stored as individual records in table 801 corresponding to that physician. The questionIDs are used to retrieve the question text

and response format from question pool table <u>808</u>. The description of the response format is then obtained from table <u>809</u>, and the question is presented to the patient. If the patient is accessing the site directly, the questions are preferably presented in the form of web pages to the patient. If data entry personnel are utilized, the question may be read to the patient to solicit the response. The responses are entered and submitted to the database system, and are stored in tables <u>802</u> and <u>803</u>. Tables <u>804</u>, <u>805</u>, <u>806</u>, <u>809</u>, <u>810</u> and <u>811</u> are used in a similar to gather compliance information.

With reference to Figure 19, tables 820, 821 and 822 are used to gather guest information, including requests for information and contact information. Tables 825 and 826 are used to provide a communication interface that facilitates communication between physicians and patients. The sender and receiver are identified, along with the subject, body, date, time and priority, etc. Thus, when a physician logs in and requests to view his email, the communication interface, preferably made up of server scripts and database management software, queries email table 826 for all entries having the physician's UserID (table 800) in either the senderID or receiverID field. Those entries may be looked up in email audit table 825 to check the status.

Table 827 keeps a log of frequently asked questions (FAQs – as stored in FAQ table 824) that have been accessed by a given patient, thus allowing a physician to obtain further information about potential concerns of his or her patients. With reference to Figure 21, Letters of consultation may be created using tables 850, 851 and 854. New diagnoses are created via the physician data access interface and stored in table 852. The template for creating a diagnosis may be created using standard diagnostic phrases obtained from tables 853 and 856. The new diagnosis template may present the cardiologist with certain diagnosis suggestions based on the test results. Specifically, the diagnostic 113 engine may perform algorithmic tests on new test results and store one or more diagnosis candidates in diagnosis fields (not shown) in the test results table 902.

With reference to Figure 22 treatment plans are stored in tables <u>860</u> - <u>865</u>. Physicians, through the physician data access interface, create treatment plans for individual patients. Typically, the treatment plan is created or updated after test results become available. The treatment plan template is preferably generated and stored in treatment plan table <u>860</u>, reading assignment table <u>861</u>, Exercise assignment table <u>862</u>,

10

15

20

25

30

diet assignment table <u>863</u>, CV product assignment table <u>864</u>, and education assignment table <u>865</u>. Tables <u>866</u>, <u>867</u>, <u>868</u>, <u>869</u>, <u>870</u>, <u>871</u> and <u>872</u> are used to store possible selections for the corresponding tables as shown in Figure 22. When a treatment plan template is generated for the physician, the diagnostic engine may be utilized to query Diagnosis table <u>852</u> and responsively determine suggested treatment plan components to be displayed as a suggested treatment plan within the template.

The diagnostic engine may also query relevant test result tables, patient history tables, compliance tables and the like to determine the suggested treatment plan. Prior to the cardiologist's review of the suggested treatment plan, the case manager or physician's assistant may initially access it to review it. The case manager may modify the plan or accept it and submit it to the cardiologist for final approval and posting to the patient's folder where the patient may view the plan.

With respect to Figure 25, the flow chart for the method 1000 will be described. At step 1002 the test result data is stored. As described herein, the test data may be received over a network from a lab. At step 1004 the test result data is supplied to the physician. Preferably, the test data is provided using web pages conveyed over a network. At step 1006 the site receives diagnostic data from the physician. The diagnostic data is stored in the database as part of the patient's record. At step 1010 the site receives treatment plan information from the physician. The treatment plan is likewise stored in the patient's record. At step 1012 the site provides treatment plan information to the patient, also via web pages conveyed over a network.

With respect to Figure 26, the flow chart for the method <u>1100</u> will be described. Steps having the same number as in the method <u>1000</u> are essentially the same. However, the method <u>1100</u> includes additional steps. At step <u>1003</u> the test result data is analyzed. The analysis is preferably performed using diagnostic engine <u>113</u>. The diagnostic analysis information is stored in a table and associated with the relevant test results. The diagnostic analysis data is preferably one or more suggested diagnoses based upon the test results and other salient data from the clinical database <u>105</u> and information database <u>104</u>. At step <u>1005</u> the test result data and the diagnostic data is supplied to the physician. At step <u>1007</u> the site analyzes the diagnostic data received from the physician. The analysis preferably identifies possible treatment plan information, including possible

prescriptions, diets, exercises, CV products, educational materials, etc., as discussed above. It should be noted that information identified as being provided to and received from the physician may actually be submitted and received by a physician's assistant, nurse, or other responsible party. The term physician is used in this context to refer to the fact that the submission of diagnosis and treatment plans should be performed under the supervision of a physician.

A preferred embodiment of the present invention has been described herein. It is to be understood, of course, that changes and modifications may be made in the embodiment without departing from the true scope of the present invention, as defined by the appended claims. The present embodiment preferably includes logic to implement the described methods in software modules as a set of computer executable software instructions. The Computer Processing Unit ("CPU") or microprocessor implements the logic that controls the operation of the site, diagnostic engine, and database. The microprocessor executes software that can be programmed by those of skill in the art to provide the described functionality.

The software can be represented as a sequence of binary bits maintained on a computer readable medium including magnetic disks, optical disks, and any other volatile or (e.g., Random Access memory ("RAM")) non-volatile firmware (e.g., Read Only Memory ("ROM")) storage system readable by the CPU. The memory locations where data bits are maintained also include physical locations that have particular electrical, magnetic, optical, or organic properties corresponding to the stored data bits. The software instructions are executed as data bits by the CPU with a memory system causing a transformation of the electrical signal representation, and the maintenance of data bits at memory locations in the memory system to thereby reconfigure or otherwise alter the unit's operation. The executable software code may implement, for example, the methods as described above.

It should be understood that the programs, processes, methods and apparatus described herein are not related or limited to any particular type of computer or network apparatus (hardware or software), unless indicated otherwise. Various types of general purpose or specialized computer apparatus may be used with or perform operations in accordance with the teachings described herein.

10

In view of the wide variety of embodiments to which the principles of the present invention can be applied, it should be understood that the illustrated embodiments are exemplary only, and should not be taken as limiting the scope of the present invention. For example, the steps of the flow diagrams may be taken in sequences other than those described, and more or fewer elements may be used in the block diagrams.

It should be understood that a hardware embodiment may take a variety of different forms. The hardware may be implemented as one or more server computers. For example, the web server may reside on one computing platform that accesses the database system residing on another computing platform.

The claims should not be read as limited to the described order of elements unless stated to that effect. In addition, use of the term "means" in any claim is intended to invoke 35 U.S.C. §112, paragraph 6, and any claim without the word "means" is not so intended. Therefore, all embodiments that come within the scope and spirit of the following claims and equivalents thereto are claimed as the invention.

ABSTRACT OF THE INVENTION

The invention relates to a cardiovascular healthcare management system which as an infomediary site with databases having information relating to cardiovascular disease risk factors such as age, blood pressure, LDL, HDL and subfractions thereof and cardiovascular disease management such as diet, exercise, drugs and cardiovascular educational materials. The physician can communicate electronically with the infomediary site to obtain patient test results and formulate a patient treatment plan from the diet, exercise and drug data. The physician can electronically communicate treatment plans to the patient through a record in the infomediary site. The patient can electronically communicate compliance information to the physician through the patient record. The infomediary site may also provide a case manager to provide initial draft treatment plans to the physician, gather patient history data and/or patient treatment plan compliance data.

1

2

3

1

2

3

2

3

4

5

6

7

8

9

1

What is claimed is:

		4 4.4			
1 4	Cordinagecular	· healthcare	management	evietem	compriging
I. <i>F</i>	A cardiovascular	meanmeare	managomon	SYSICILL	comprising.

- (a) an infomediary site having databases for cardiovascular healthcare management;
- (b) a data entry interface for receiving patient personal data and test results and storing the data and results in the infomediary site databases;
 - (c) a diagnostic engine for analyzing patient test results;
- (d) a physician data access interface to allow physician access to the infomediary databases; and
- (e) a communication system allowing the physician to communicate cardiovascular healthcare management information to the patient.
- 2. The cardiovascular healthcare management system of claim 1 further comprising a cardiovascular knowledge base that stores information related to cardiovascular risk factors.
- 3. The cardiovascular healthcare management system of claim 1 wherein the diagnostic engine includes algorithms for associating test results with possible treatments.
- 4. The cardiovascular healthcare management system of claim 1 wherein the diagnostic engine includes algorithms for associating test results with possible diagnoses.
- 5. The cardiovascular healthcare management system of claim wherein the diagnostic engine includes algorithms for associating diagnosis information with possible treatment plans.
- 6. The cardiovascular healthcare management system of claim 1 wherein the diagnostic engine for comparing patient's test results compares inputted patient LDL subfraction and HDL subfraction values to normal LDL subfraction and HDL subfraction values.
- 7. The cardiovascular healthcare management system of claim 6 wherein the normal LDL subfraction and HDL subfraction values are determined base in part of the values stored in the infomediary database.

2

1

2

3

- 1 8. The cardiovascular healthcare management system of claim 6 wherein the 2 physician dynamically selects parameters for patient test results comparisons.
- 1 9. The cardiovascular healthcare management system of claim 1, wherein the patient provides compliance data that is stored in the records for later review by the physician.
 - 10. The cardiovascular healthcare management system of claim 1, further comprising a patient access interface whereby the patient accesses a cardiovascular treatment plan.
- 1 11. The cardiovascular healthcare management system of claim 1, wherein the diagnostic engine analyzes test results and provides suggested diagnoses to the physician.
 - 12. The cardiovascular healthcare management system of claim 1, wherein the diagnostic engine analyzes diagnostic information and provides suggested treatment plans.
 - 13. A method of managing cardiovascular healthcare information, comprising the steps of:
 - (a) storing test result data in a database;
 - (b) providing test result data to a physician via a network;
 - (c) receiving diagnostic information from a physician via a network;
 - (d) receiving treatment plan information from a physician via a network; and
 - (e) providing treatment plan information to a patient via a network.
- 1 14. The method of managing cardiovascular healthcare management system of claim 2 13, wherein the step of providing test result data is performed using internet protocols.
- 1 15. The method of managing cardiovascular healthcare management system of claim 2 13, wherein the step of receiving diagnostic information from a physician is performed using 3 internet protocols.
- 1 16. The method of managing cardiovascular healthcare management system of claim 2 13, wherein all the steps of receiving and providing information are performed using internet 3 protocols.

- 1 17. The method of managing cardiovascular healthcare management system of claim
- 2 13, further comprising the steps of analyzing the test result data and providing diagnosis
- 3 information to the physician.
- 1 18. The method of managing cardiovascular healthcare management system of claim
- 2 13, further comprising the steps of analyzing the received diagnosis information and providing
- 3 suggested treatment plan information to the physician.
- 1 19. The method of managing cardiovascular healthcare management system of claim
- 2 13, further comprising the step of receiving compliance information from a patient.
- 1 20. The method of managing cardiovascular healthcare management system of claim
 - 13, further comprising the step of receiving compliance information from a patient.





FIG. 1



Submit Guest

150













Figure 4






Contact L

About 4myheart.com





Physicians Front Desk

Physicians Front Desk

About 4myheart.com

Contact

personalized cardiovascular wellness management

Patient Folder

WELCOME
Patient Activity
Patient List
E-Mail
My Statistics
End of Day Report
Preferences
Reference
News
Product Review

Patient Joeseph Patient

Gender Joeseph Faria

Age 57
Birth Date 4-16-4

Birth Date 4-16-42 ID Number T25048

Salient Medical History Current Medications

Information Information

Test Date 05-28-1999

HDL Gradient Gel Electrophoresis (HDL GGE); Relative HDL Subclass Distribution

HDL Subclass
2b
2a
Sum of HDL2

results (area %)	10%	50%	90%
15	4	8	18
22	12	18	26
37	Sound of the Carte		
(and what has been been been been the second of the contract o	SANDONE OF THE PROPERTY OF THE		Роментинго и постано по выполняющих постановах по выполняющих постановах по выполняющих постановах по выполняющих по

HDL Subclass
350000000000000000000000000000000000000
3a
Bayer-manufacture production of the contract o
3b
Engineering continue and a single and an analyze and a supplementary and a supplementa
3c
Sum of HDL3
CONTRACTOR OF THE PROPERTY OF

-	results (area %)	10%	50%	90%
0.151001010101	32	31	36	40
The state of the state of	23	16	26	33
	9	5	10	17
A STATE OF THE PERSON NAMED IN	64			

HDL2 and HDL3 contribute to HDL-cholestetol measurements. HDL3 levels are more resistant to change than HDL2, which can be increased by weightloss, excercise and estrogen.

Particle Size (numbers at peaks are particle size in angstroms)

HH0521.BIO:8 B99006263 2G






Figure 8




The state of the s

260



Contact U





Physicians Front Desk About 4myheart.com

· personalized cardiovascular wellness management

WELCOME
Patient Activity
Patient List
E-Mail
My Statistics
End of Day Report
Preferences
3644 CEEE
Reference
News
Product Review
A hand chart minimum as a second second

ruiovascular welliness	The state of the s				
Patient Folder	1000 pt		(Dia	gnosis Treatmer	t Plan Letter of Cons
Patient	Joeseph Patio	ent			
Gender	Male				
Age	57				1 4 5
Birth Date	4-16-42				
ID Number	T25048				
Salient Medical Hi	story	C	Current Medication	ons	
Information		i	nformation		
Overview 30	2				320 -Q
· · · · · · · · · · · · · · · · · · ·	paragraph explaining t	he comments	of the treatment p	lan below.	11
Exercise - 304	1			**************************************	\\\a
High Caloric Burn					
Prescription Medic	-306	utadouandodeleMilliodorenhimentos			T da
Drug 1 - 100mg - Q	AM - 250 count - 2 Re	efills - No Food			
Diet - 308			***************************************	nuu annuu muun makkii kali, Mada palaadi akkii	<u> </u>
3000 Total Calories	;				
15% Polyun	ted Fat saturated Fat insaturated Fat		10% 45% 15%	Complex Carboh Simple Carbohyo Protiens	
CV Products - 3	310				
Product					
Product					
Product Product					
Suggested Readin	310	T			- La
Title	igs - JIZ	Туре			
Title		Book			
Title	,	Journal			
http://www.website.	.com	Website			
Adjunct Medication	-314				- Q
Name					
Follow Up	316				\ a
Next Office Visit		Next Sche	duled Test(s)		
3 Months			LDL, Lipid Pane	I	£





Physicians Front Desk

About 4myheart.com

Contact L

personalized cardiovascular wellness management



Patient Folder		Diagnosis Treatment Plan Letter of Cons
Patient	Joeseph Patient	
Gender	Male	
Age	57	
Birth Date	4-16-42	
ID Number	T25048	

Salient Medical History Current Medications

Information Information

Diet Prescription

Total Daily Calories


Select Daily Caloric Amount 2000 3000 4000


Total Caloric Breakdown


6% 15% 10%	Saturated Fat Polyunsaturated Fat Monounsaturated Fat	10% 45% 15%	Complex Carbohydrates Simple Carbohydrates Protiens
6%	Saturated Fat	10%	Complex Carbohydrates
15%	Polyunsaturated Fat	45%	Simple Carbohydrates
10%	Monounsaturated Fat	15%	Protiens
6%	Saturated Fat	10%	Complex Carbohydrates
15%	Polyunsaturated Fat	45%	Simple Carbohydrates
10%	Monounsaturated Fat	15%	Protiens
6%	Saturated Fat	10%	Complex Carbohydrates
15%	Polyunsaturated Fat	45%	Simple Carbohydrates
10%	Monounsaturated Fat	15%	Protiens

Notes









Patient Reception Area

About 4myheart.com

Contact L

personalized cardiovascular wellness management



H. Robert Superko

Requesting Physician

Research Study - Lipid Institute 1875 South Grant Street, Suite 700 San Mateo, CA 94402 FAX (650) 372-1952

Patient

D*****, Frank Gender

Age

ID Number T14523

Birth Date 10-20-1927

Report Details

Test Date

05-28-1999

HDL Gradient Gel Electrophoresis (HDL GGE) : Relative HDL Subclass Distribution

HDL Subclass	200000000000000000000000000000000000000
2b	
2a	7
Sum of HDL2	


57.7				
4		garanta aranda transcribita de la companya de la c	process and the process of the proce	Contraction recommendation of the contraction of th
3	ensulta (area 21)	109	EUው	፤ ዐ∩ጭ ∦
H	results (area %)	1040	30%	30%
- 8		janen er en er		
9	10		0	10 1
8	10 1	4		10 1
200	CONTRACTOR OF THE PROPERTY OF		NORMAN CONTRACTOR AND A	CONTRACTOR OF THE PROPERTY OF
- 8	77	17	10	70 1
200		1 14 1	10 1	_ ZD 1
-	AND THE RESIDENCE OF THE PARTY			
-	ኅማ			is in
- 5	J/ 1			9


HDL Subclass
LINE AMERICA
Experimental control of the control
3a
1 JO
President and the second and the sec
3b
i 3D
parameter and the second secon
1 7_
3c
TO COMPAN
Sum of HDL3

results (area %)	10%	50%	90%
32	31	36	40
23	16	26	33
9	5	10	17
64			

HDL2 and HDL3 contribute to HDL-cholestetol measurements. HDL3 levels are more resistant to change than HDL2, which can be increased by weightloss, excercise and estrogen.

Particle Size (numbers at peaks are particle size in angstroms)









Figure 16






Figure 18

tGuestInfoResponseItemis ∭ ID

RequestID ResponseTypeID

Figure 20

NORTHGUE NEW TOUR






Figure 21




Cardiovascular Informatics

				Fibrinogen ACAT PPARA	CETP LPL MTHFR	Apo E Apo B
		` ;	900	MTP FDFT1 HMGCR	CBS GPIII	PON ACE AGT
Niacir Niacir Niacir Folate Folate Hesin Low ft				902 -	1	
POSS Niacin-Fibrates-Diet Niacin-Estrogen Folate-Pyridoxine (B6) Resin-Statin-Niacin Low fat diel-Probucol + Niacin-Fibrates		Car Therap that sy		Apo B Apo E Lp(a)	Homocysteine Fibrinogen Insulin	Total Cholesterol Total HDL Total LDL Triglycerides
POSSIBLE TREATMENTS s-Diet Fenofibrates en MTP inhibitors Estrog ine (B6) PPARA inhibitors Do no Viacin ACAT inhibitors Vitami vibrates Niacin-Fibrates ALEAN TO THE TRANSPORT OF THE TRANSPORT	(=)	Cardiovascular Therapeutic Database that synthesizes data	(#)		TDLQ12 TDLQ18	TDFO14 TDFO12 TDFO11
MENTS as itors hibitors hibitors hibitors rates		base data		Peak perlicle size (dm)	HDLOTS IDLOTS	HDLQ71 HDLQ73 HDLQ73
Resins Estrogen Do nothing Vitamin E Fibrates		906		704	1	•
908		Hawre	Hypothyroidsm Renal failure	Hypertholesterohemia Hypertipidemia Hyperthyroidism hypertivoemia	Obesity Peripheral vescular disease Homocysteinemia	Angina Coronary aftery disease Stroke Hypertension
	((ó			50	

UDIECO" WITHERD





Case No.: MBHB00-203

DECLARATION AND POWER OF ATTORNEY FOR PATENT APPLICATION

As a below named inventor, I hereby declare that:

My residence, post office address and citizenship are as stated below next to my name.

I believe I am the original, first and sole inventor (if only one name is listed below) or an original, first and joint inventor (if plural names are listed below) of the subject matter which is claimed and for which a patent is sought on the invention entitled:

C	CARDIOVASCULAR I	HEALTHCARE MANA	GEMENT SYSTEM AND METHOD	
the specification	on of which is attached l	nereto unless the followin	ng space is checked:	
	was filed on	as United	States Application Serial Number	
I hereby state claims, as ame	that I have reviewed an	d understand the content treferred to above.	ts of the above-identified specification, inclu-	ding the
I acknowledge	e the duty to disclose inf	ormation which is mater	ial to patentability as defined in 37 CFR § 1.50	5.
patent or invector country other application for	entor's certificate, or § than the United States, !	365(a) of any PCT into disted below and have also tificate, or PCT internation	9(a)-(d) or § 365(b) of any foreign application ernational application which designated at less identified below, by checking the box, any onal application having a filing date before the	east one foreign
Prior Foreign Num 1.	Application(s): ber Country	Day/Mor	nth/Year Filed	
<u>Appl</u>	n the benefit under 35 U. lication Number 126,003	S.C. § 119(e) of any Un Filing Date March 24, 199	ited States provisional application(s) listed bel	.ow:
	168,354	December 1, 1		
international a claims of this provided by the to patentabilia application and	application designating to application is not discle the first paragraph of 35 ty as defined in 37 C.	the United States, listed losed in the prior United U.S.C. § 112, I acknowle	United States application(s), or § 365(c) of a pelow and, insofar as the subject matter of each States or PCT international application in the edge the duty to disclose information which is same available between the filing date of this application. Status: patented, pending, abandone	ch of the manner material the prior
	oint the following attorn this application and		business in the Patent and Trademark connected therewith:	k Office

Full name of first inventor:

Denis A. Berntsen	Reg. No. 26707	George I. Lee	Reg. No. 39269
John J. McDonnell	Reg. No. 26949	James M McCarthy	Reg. No. 39296
Daniel A. Boehnen	Reg. No. 28399	Monica Grewal	Reg. No. 40056
Bradley J. Hulbert	Reg. No. 30130	Jeremy Noe (agent)	Reg. No. 40104
•	Reg. No. 30243	Sean M. Sullivan	Reg. No. 40191
Paul H. Berghoff	_	Timothy R. Baumann	Reg. No. 40502
Grantland G. Drutchas	Reg. No. 32565	•	U
Steven J. Sarussi	Reg. No. 32784	Amir N. Penn	Reg. No. 40767
David M. Frischkorn	Reg. No. 32833	Patrick J. Halloran	Reg. No. 41053
James C. Gumina	Reg. No. 32898	Joshua R. Rich	Reg. No. 41269
A. Blair Hughes	Reg. No. 32901	Thomas E. Wettermann	Reg. No. 41523
Thomas A. Fairhall	Reg. No. 34591	Vernon W. Francissen	Reg. No. 41762
Emily Miao	Reg. No. 35285	Robert J. Irvine	Reg. No. 41865
Kevin E. Noonan	Reg. No. 35303	Richard A. Machonkin	Reg. No. 41962
Leif R. Sigmond, Jr.	Reg. No. 35680	David S. Harper	Reg. No. 42636
Lawrence H. Aaronson	Reg. No. 35818	Stephen Lesavich	Reg. No. 43749
Matthew J. Sampson	Reg. No. 35999	Enrique Perez	Reg. No. 43853
Curt J. Whitenack	Reg. No. 36054	Marcus J. Thymian	Reg. No. 43954
Christopher M. Cavan	Reg. No. 36475	S. Richard Carden	Reg. No. 44588
Michael S. Greenfield	Reg. No. 37142	Mark L. Chael (agent)	Reg. No. 44601
Roger P. Zimmerman	Reg. No. 38670	Stephen H. Docter	Reg. No. 44659
Anthoula Pomrening	Reg. No. 38805		

Address all telephone calls to	Robert J. Irvine III	at (312) 913-0001.

Frank R. Ruderman

Post Office Address: 419 Baden Street, San Francisco, California 94131-2831

Address all correspondence to McDonnell Boehnen Hulbert & Berghoff, 300 South Wacker Drive, Chicago, Illinois 60606 USA.

I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code and that such willful false statements may jeopardize the validity of the application or any patent issued thereon.

Inventor's signature:	Date:
Residence: 248 Bayview Drive	
San Carlos, California 94070	
Citizenship: United States	
Post Office Address: 248 Bayview Drive	
San Carlos, California 94070	
Full name of second joint inventor: David T. Shewmake	
Inventor's signature:	Date:
Residence: 419 Baden Street, San Francisco, California 94131-2831	
Citizenship, United States	